	Product Development	Clause Ref : ISO 13485:2003;ISO 9001:2008
		Version : 1.0
	InstaRISPACS-IHE Conformance Statements	Effective Date : 05/03/2013
		Doc # : InstaRISPACS-IHE Integration Statement

InstaRISPACS IHE Integration Statement

	Product Development	Clause Ref : ISO 13485:2003;ISO 9001:2008
		Version : 1.0
	InstaRISPACS-IHE Conformance Statements	Effective Date : 05/03/2013
		Doc # : InstaRISPACS-IHE Integration Statement

Reference Documents:


Sl. No.	Doc #

Revision History:

SI No.	Revision#	Date	Description	Author	Approver
1	1.0	1-Mar-2012	Version 1.0	Angad Nath	Tapesh Kumar Agarwal
2	1.1	9 th July 2014	Version 1.1	Tapesh K Agarwal	
3	1.2	16 th May 2016	Removing 'Report Creator' from Reporting Workflow (RWF) profile	Tapesh K Agarwal	
4	1.3	10 th Sept 2017	Adding the profile IOCM	Tapesh K Agarwal	

Disclaimer:

This document is intended only for the person to whom it is addressed and/or otherwise authorized personnel. The information contained herein and attached is confidential and the property of Meddiff Technologies Pvt. Ltd.


	Product Development	Clause Ref : ISO 13485:2003;ISO 9001:2008
		Version : 1.0
	InstaRIS PACS-IHE Conformance Statements	Effective Date : 05/03/2013
		Doc # : InstaRIS PACS-IHE Integration Statement

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profile	Actors (Implemented)	Options (Implemented)
Access to Radiology Information (ARI)	1. Image Manager 2. Image Archive 3. Image Display	None None None
Scheduled Workflow (SWF)	1. Image Manager 2. Image Archive 3. Performed Procedure Step Manager 4. Image Display 5. Order Filler 6. Order Placer	None None None None None None
Consistent Presentation of Images (CPI)	1. Image Manager 2. Image Archive 3. Image Display	None None None
Patient Information Reconciliation (PIR)	1. Image Manager 2. Image Archive 3. Performed Procedure Step Manager	None None None
Portable data for Imaging (PDI)	1. Image Display 2. Portable Media Importer 3. Portable Media Creator 4. Report Reader	None None None None
Cross Enterprise Document Sharing for Imaging	1. Document Registry 2. Document Repository 3. Document Consumer 4. Imaging Document Source 5. Imaging Document Consumer	None None None Set Of DICOM Instances and Text Report None
Audit Trail and Node Authentication (ATNA)	1. Secure Application	None
Presentation of Grouped Procedures (PGP)	1. Image Archive 2. Image Manager	None
Key Image Note (KIN)	1. Evidence Creator 2. Image Display 3. Image Manager 4. Image Archive	None
Post-Processing Workflow (PWF)	1. Image Display 2. Image Manager 3. Image Archive	None

	Product Development	Clause Ref : ISO 13485:2003;ISO 9001:2008
		Version : 1.0
	InstaRISPACS-IHE Conformance Statements	Effective Date : 05/03/2013
		Doc # : InstaRISPACS-IHE Integration Statement

	4.Evidence Creator	
Reporting Workflow (RWF)	1.Report Manager 2.Report Creator 3.Image Manager 4.Image Archive	None
Evidence Document (ED)	1.Image Manager 2.Image Archive	None
Basic Security (SEC)	<i>Deprecated</i>	
Consistent Time Client (CT)	1.Time Client	None
Cardiac Catheterization Workflow (CATH)	1.Image Display 2.Image Manager 3.Image Archive	None
Simple Image & Numeric Report (SINR)	1.Report Creator 2.Report Manager 3.Report Repository 4.Report Reader	None
Radiation Exposure Monitoring (REM)	1.Image Archive 2.Image Manager 3.Dose Information Consumer	None
Mammography Image (MAMMO)	1. Image Manager 2. Image Archive 3. Image Display	None
Charge Posting (CHG)	1. Order Filler 2. ADT Patient Registration	None
PIX – Patient Identifier Cross-referencing	Patient Identifier Cross-reference Consumer	None
XDS.b – Cross Enterprise Document Sharing	1. Document Repository 2. Document Consumer	None
XDS-I.b – Cross Enterprise Document Sharing for Imaging	1. Document Repository 2. Imaging Document Source 3. Imaging Document Consumer	None
PIR – Patient Information Reconciliation	1. Image Manager / Archive 2. Performed Procedure Step Manager	None
PDQ – Patient Demographics Query	1. Patient Demographics Consumer	None
XDS-SD – Cross Enterprise Sharing Of Scanned Documents	1. Document Consumer	None
Image Object Change Management (IOCM)	1. Image Manager 2. Image Archive 3. Change Requestor	None

	Product Development	Clause Ref : ISO 13485:2003;ISO 9001:2008
		Version : 1.0
	InstaRISPACS-IHE Conformance Statements	Effective Date : 05/03/2013
		Doc # : InstaRISPACS-IHE Integration Statement